EASY SMILE

PRO

MANUAL DE USO

KIT DE BLANQUEAMIENTO DENTAL

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Certificate of Registration

This is to certify that the Quality Management System of

Nanchang Dental Right Technology Co., Ltd.

Unifast Social Credit Code: 913601060702048P

Operation Address: Building 101, Building 102, Nanchang Eastern Industrial Park, No.2799,
Tiansheng Road, Hupu-Sao Zone, Nanchang City, Jiangxi Province, China

Registered Address: Building 101, Nanchang Jiefang Industrial Park, No.2799, Tiansheng Road,
Hupu-Tech Zone, Nanchang City, Jiangxi Province, China

applicable to

Production and sales of tooth whitening instruments, tooth whitening gels and tooth whitening kit (export to Europe and America)

has been assessed and approved by NQA against the provisions of

ISO 13485:2016

This organization is subject to be monitored according a quality management system, to the above standard, which will be awarded by NQA.

Certified Client shall accept regular surveillance assessments, the validity of certificate shall be maintained for the positive result of audit.

The certification of this certificate can be checked on CNCA's website (www.cnca.gov.cn)

NQA's website: www.nqa.com.cn

Managing Director

Certificate Number: 47344

Date: 22 December 2019
Valid Until: 22 December 2024
EAC Code: 18

CE

Signed for and on behalf of
 ILL-CRSC Standards Technical Services Co., Ltd. (Singapore Branch)

Jenny
Jenny Lee
Approved Signature

Downloaded from ascelibrary.org by University of California, San Diego on 06/01/14. Copyright ASCE, For All Rights Reserved, No part of this document may be reproduced without written permission from ASCE.

[illegible]

	
Verification Report	
No. CASBC/CH220101	Date: 18 Aug 2016
Page 1 of 23	
BUCHANG TECH. BRIGHT TECHNOLOGY CO., LTD. BUCHANG 2ND. DACHANGJIAN ZHANG INDUSTRIAL PARK, 279 THUNDERING HWY., HONG-TECH ZONE, SUICHANG, JIANGSU, CHINA	
Sample Name :	Tenth Wharfing Lead Light
SGS Job No. :	CPN-04700 - 002
Internal Reference No. :	CPN-04700001
Tested Items Under Test :	L1-L102
Test Method / Standard :	ISO 9001
Manufacturer :	L1-L107
Manufacturer :	BUCHANG TECH. BRIGHT TECHNOLOGY CO., LTD.
Date of Sample Received :	25 JUN 2016
Verification Period :	25 JUN 2016 - 18 AUG 2016
Verification Requested :	Audit according to Quality Certificate (20140805)-renewing ISO9001
Verification Method :	Please refer to the audit plan.
Verification Report :	Please refer to the audit plan.
Verification Conclusion :	Based on the verification results of the submitted samples, the results of audit, observation, consultation, discussion, Pre-qualified Inspection Report (PPQR), Pre-qualified Assessment (APR) together with the limits set and by BUCHANG 2ND. DACHANGJIAN ZHANG INDUSTRIAL PARK, 279 THUNDERING HWY., HONG-TECH ZONE, SUICHANG, JIANGSU, CHINA.
Note :	The test results are related only to the tested items. The report shall not be considered as evidence in full without the entire assessment of the factory facilities.

Signal for and on behalf of
S&S/CSTC Standards Technical Services Co., Ltd. Quartermen Road

Jenny Liao
Approved Signature

© 2006 Blackwell Publishing Ltd, *Journal of Internal Medicine* 260: 105–114

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